



PHILIPS

IntelliSite

Pathology Solution

Remote use

Facilitate digital primary diagnosis via expansion of remote use*

Specifications USA

FDA's enforcement discretion allows for the use of readily available consumer monitors at home with the Philips IntelliSite Pathology Solution during COVID-19 emergency.

The solution is designed for in-vitro diagnostic purposes, helping pathologists to review and interpret digital images of surgical pathology slides prepared from formalin-fixed paraffin-embedded (FFPE) tissue samples. It ultimately helps your team improve patient care by enhancing knowledge-sharing and maximizing resources

More productivity through digital case sign-out is possible today with the Philips IntelliSite Pathology Solution for routine histopathology diagnosis. We can help you throughout your digital pathology transformation journey with the aim to improve your decision-making in precision diagnostics.

Philips IntelliSite Pathology Solution (PIPS) for remote use

The specified device and use are according to FDA's enforcement discretion allowance

Scanner	Ultra Fast Scanner (UFS)	
Access to image server	Through Virtual Private Network (VPN)	
Viewing software	IntelliSite Web-based browser	
Client Computer Specifications	CPU	Intel processors with more than two (2) cores and 2.3GHz (Intel i5 or above is preferred)
	RAM	3GB (8GB or above is preferred) Hard Drive – 50GB (250GB or above is preferred)
	Operating System	Windows 7/10 (Win 10 Pro is preferred)
	Video monitor connection	Graphical card with connection for two monitors is preferred
Internet connection speed	100Mbps or above	
Monitor (Display)	Brightness	≥ 300 Cd/m2;
	Contrast Ratio	>1000:1
	Color Gamut	>99% sRGB
	Color temperature	D65 or similar
	Delta E	<3 or similar
	Display technology	IPS with LED backlight
		Recommended models: Philips PP27QHD (MMPC-4227F1) (FDA Cleared PIPS 3.2) Philips PS27QHDCR (MMCC-4127) Barco MDRC-2324 Barco Nio MDNC-3421 Eizo RX240 Eizo MX216 Eizo MX242W Eizo CS2731 Dell MR2415 Dell MR2416 Dell UP3017 Dell U2415 Dell UP2716D Dell UP2719

*The expansion of remote use of Philips IntelliSite Pathology Solution in the USA are under specific conditions and device specifications. The type of monitors is only allowed due to emergency situations. The limitations and warnings specified in device labeling are still applicable to this situation.

This material is only for distribution in the USA

© 2020 Koninklijke Philips N.V. All rights reserved.
 Specifications are subject to change without notice.
 Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



www.philips.com/digitalpathology

4522 991 58981 * APR 2020