





|                               |                                   |  |
|-------------------------------|-----------------------------------|--|
| <b>PHILIPS<br/>HEALTHCARE</b> | <b>Product Recycling Passport</b> |  |
|-------------------------------|-----------------------------------|--|


|                                   |                             |   |
|-----------------------------------|-----------------------------|---|
| <b>Product name</b>               | <b>Volcano s5 PC Module</b> |   |
| <b>Identification code(s)</b>     | <b>807300-001</b>           |   |
| <b>Total weight (in Kg)</b>       | <b>70</b>                   |   |
| <b>Producer/<br/>Manufacturer</b> | <b>Name<br/>company:</b>    | <b>Philips Volcano</b>                      |
|                                   | <b>Address:</b>             | <b>2870 Kilgore Road Rancho Cordova, CA</b> |
|                                   | <b>Zip code:</b>            | <b>95670</b>                                |
|                                   | <b>Country:</b>             | <b>USA</b>                                  |


| <b>Recycle Info</b>  | <b>Items:</b> | <b>Location</b> |
|--|---------------|-----------------|
| <b>Special attention</b><br><br> | <b>N/A</b>    |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |

|  | <b>Items:</b> | <b>Location</b> |
|--|---------------|-----------------|
| <br><br><b>Fluids / Gases</b> | <b>N/A</b>    |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |
|  |               |                 |

|   | <b>Type:</b>                                   | <b>Location</b>    |
|---|--|--------------------|
| <br><br><b>Batteries<br/>To be Removed</b> | <b>CR2032 Button Cell Lithium Battery 3.0V</b> | <b>Motherboard</b> |
|   |  |                    |
|   |  |                    |
|   |  |                    |
|   |  |                    |

|   | Substances:                     | Location          |
|---|---------------------------------|-------------------|
| Hazardous   |                                 |                   |
|  | may contain lead                | Circuit board     |
|   | may contain hexavalent chromium | Screws, fasteners |
|   |                                 |                   |
|   |                                 |                   |
|   |                                 |                   |
| To be Removed   |                                 |                   |





Motherboard Location



Motherboard Battery to be removed



### Recycling Passport Revision History

| Rev | Description       | Date   |
|-----|-------------------|--------|
| A   | S5x WEEE Passport | 4/7/16 |