





PHILIPS HEALTHCARE	Product Recycling Passport	
-------------------------------	-----------------------------------	--


Product name	Volcano SyncVision System	
Identification code(s)	400-0100.10	
Total weight (in Kg)	70	
Producer/ Manufacturer	Name company:	Philips Volcano
	Address:	2870 Kilgore Road Rancho Cordova, CA
	Zip code:	95670
	Country:	USA

Recycle Info	Items:	Location
Special attention 		
	N/A	

 Fluids / Gases		

 Batteries To be Removed		

Hazardous	Substances:	Location
	may contain lead	Circuit board
	may contain hexavalent chromium	Screws, fasteners
To be Removed		





**Coin type motherboard
battery to be removed**



Recycling Passport Revision History

Rev	Description	Date
A	SyncVision WEEE Passport	4/7/16